

TONGUE & LIP PRESS



**HARD VERSION
TLP-H**




myoTALIA[®]

**PATENTS PENDING WORLDWIDE,
FOR MORE INFORMATION VISIT MYORESEARCH.COM**

INDICATIONS FOR USE

The *Myotalea*[®] TLP is to be used as an aid in alleviating symptoms. A simple sequence will assist in developing better muscle control and build up to the required time durations based on individual tolerance to fatigue in order to build strength and tone.



TONGUE

Tongue Press & Hold

- Position the press tube between the tongue and the roof of the mouth.
- Compress the entire body of the tongue firmly against the press tube into the roof of the mouth. Hold for 5 breaths while keeping lips together and breathing lightly through the nose.
- See below for an advanced technique.*


Active Tongue Press

- Position the press tube between the tongue and the roof of the mouth. The teeth can be slightly separated while doing this exercise.
- While maintaining a light nasal breathing pattern, elevate the entire body of your tongue and actively compress the press tube closed into the roof of your mouth for 1 breath, and release. Repeat this step 5 times.
- See below for an advanced technique.*

***HEAD TILT VARIATION:** An advanced technique can be used for an interval of 5 breaths to further stretch and tone.

Use the softer TLP-S if you are having difficulty.

ating sleep-related breathing disorders; the muscle strength and tone. It is recommended to current functional abilities and push the muscles



LIPS

Lip Press & Hold

- Position the middle of the *TLP* between the lips, just in front of the teeth.
- Maintain tongue in the roof of the mouth and press lips firmly together to compress the press tube closed. Hold for 5 breaths while breathing lightly through the nose.
- See below for an advanced technique.*

Active Lip Press

- Position the middle of the *TLP* between the lips, just in front of the teeth. The teeth can be slightly separated while doing this exercise.
- Maintain tongue in the roof of the mouth while breathing lightly through the nose. Press lips firmly together to actively compress the press tube closed and hold for 1 breath, then release. Repeat this step 5 times.
- See below for an advanced technique.*

performed by gently tilting the head back and holding get the oropharyngeal muscles of the head and neck.

ultly completing exercises with the *TLP-H*.

DESIGN FEATURES


- Handle with strap tether hole can be used to attach to other devices.
- Strengthens the tongue and lips, allowing for improved muscle function.
- Made from medical grade silicone and available in two strength levels.

Press tube



Flow Sequence

If advised by your practitioner, these exercises can be combined to form an exercise flow sequence as follows:



USE FOR
5 MINUTES
TWICE
PER DAY



Handle

Strap
tether
hole



TONGUE & LIP PRESS

Head Tilt Variation

Repeat below sequence with **head tilt posture**. There should be no discomfort in head, neck or jaws during this variation.

Lip Press
& Hold
(5 breaths)




Active
Lip Press
(5 sets)



Lip Press
& Hold
(5 breaths)

The **Myotalea® Tongue & Lip Press (TLP)** has been designed to strengthen and tone the key muscles of the mouth and lips. Furthermore, the *Myotalea® TLP* targets the oropharyngeal muscles that are directly involved in keeping the upper airway open, working to improve its patency during sleep.

The *Myotalea® TLP* can be used with a head tilt postural variation to directly target and engage the key muscles of the neck and throat that are prone to collapsing during sleep-related breathing disorders.



Tongue Press

myoTALEA®
TONGUE AND LIP
EXERCISE APPLIANCE

Myotalea[®] is a range of myofunctional appliances that are used to assist in the correction of poor oral habits. The appliances can be used in conjunction with the *Myobrace*[®] and *Myosa*[®] treatment systems. The *Myotalea*[®] TLP helps to increase muscle strength, maintain the correct tongue resting position and establish adequate lip seal. This appliance has two firm levels - soft and hard, allowing for a smoother progression as muscle strength develops. Exercises are used to achieve the required outcome.

The *Myotalea*[®] TLP is designed specifically to improve lip and tongue strength, which aids in the correction of poor oral habits, including:

- Mouth breathing
- Low tongue resting position
- Incorrect swallowing pattern
- Lips apart at rest



Correct breathing



Correct swallowing

CLEANING INSTRUCTIONS: Clean your *Myotalea*[®] TLP appliance under warm running water every time you remove it from your mouth. Use *Myoclean*[™] tablets to correctly clean twice every week. Ask your doctor or dental therapist about *Myoclean*[™], the recommended cleaning agent for all *MRC* appliances.

MYOTALEA® TLP

Myotalea® Tongue & Lip Press (TLP) by *Myofunctional Research Co. (MRC)* is an active myofunctional appliance, which targets the key muscles of the tongue and lips that may require additional strength and tone in order to obtain a stable orthodontic result. Poor muscle strength/tone in the tongue and lips can also be a contributing factor to increased risks for sleep-related breathing disorders.

CAUTION: CHILDREN MUST BE SUPERVISED BY AN ADULT WHEN USING THE MYOTALEA® TLP.

Lifetime use: 6 months based on everyday wear.
Replace after 12 months regardless of wear time.

HARD version (TLP-H)



**MYOFUNCTIONAL
RESEARCH CO.**

INNOVATIVE DENTAL APPLIANCE
TECHNOLOGY AND EDUCATION



Myotalea® by *Myofunctional Research Co.*
Myofunctional Research Co. Australia
44 Siganto Drive, Helensvale QLD 4212
Australia Tel: +61 7 5573 5999

USA Representative:

Myofunctional Research Co. USA
9267 Charles Smith Avenue,
Rancho Cucamonga CA 91730 USA
PO Box 2781, Rancho Cucamonga CA 91730
USA Tel: +1 866 550 4696



Europe Representative:

Myofunctional Research Europe B.V.
Gompenstraat 21c 5145 RM Waalwijk
P.O. Box 718, 5140 AS Waalwijk
The Netherlands Tel: +31 416 651 696

Single Patient Use.

