

Specifically designed for the
Mouth Breather and Snorer



TMJ.MBV™
MOUTH BREATHER VERSION

www.myoresearch.com


WORLDWIDE PATENTS, FOR MORE INFORMATION VISIT MYORESEARCH.COM

BIGGER SIZE FOR BETTER RETENTION

The *TMJ-MBV* is specifically designed for mouth breathers who suffer from TMJ disorder symptoms. The higher vertical sides of the *TMJ-MBV* provide better comfort and retention, particularly when the mouth is open.

TMJ-MBV FEATURES

- ① **Four Large Breathing Holes**
allow for easy breathing with the mouth open.
- ② **Tongue Guard**
controls the tongue position.
- ③ **Higher Sides**
for optimum retention.
- ④ **Air Spring Base**
for variable bite opening and gentle TM Joint decompression.



Cleaning the TMJ-MBV: Clean your *TMJ-MBV* under warm running water every time you remove it from your mouth. Use *Myoclean*™ tablets to correctly clean twice every week. Ask your Doctor about *Myoclean*™, the recommended cleaning agent for all *MRC* Appliances.

CAUTION: USE ONLY UNDER THE SUPERVISION OF A LICENSED ORTHODONTIST OR DENTIST

Lifetime use: 6 months based on everyday wear.
Replace after 12 months regardless of wear time.



by Myofunctional Research Co.
Myofunctional Research Co. Australia
44 Siganto Drive, Helensvale QLD 4212 Australia
Tel: +61 7 5573 5999

USA Representative:

Myofunctional Research Co. USA
9267 Charles Smith Avenue, Rancho Cucamonga CA 91730 USA
PO Box 2781, Rancho Cucamonga CA 91730 USA
Tel: +1 866 550 4696



Europe Representative:

Myofunctional Research Europe B.V.
Gompensstraat 21c 5145 RM Waalwijk The Netherlands
P.O. Box 718, 5140 AS Waalwijk The Netherlands
Tel: +31 416 651 696



**MYOFUNCTIONAL
RESEARCH CO.**

INNOVATIVE DENTAL APPLIANCE
TECHNOLOGY AND EDUCATION

Single Patient Use.



KEEP OUT
OF DIRECT
SUNLIGHT



CLEAN
BEFORE USE